

April 2020

The self-study lesson on this central service topic was developed by 3M Health Care. The lessons are administered by Endeavor Healthcare Media.

**Earn CEUs**

After careful study of the lesson, complete the examination at the end of this section. Mail the completed test and scoring fee to *Healthcare Purchasing News* for grading. We will notify you if you have a passing score of 70 percent or higher, and you will receive a certificate of completion within 30 days. Previous lessons are available at [www.hponline.com](http://www.hponline.com).

**Certification**

The CBSPD (Certification Board for Sterile Processing and Distribution) has pre-approved this in-service for one (1) contact hour for a period of five (5) years from the date of original publication. Successful completion of the lesson and post test must be documented by facility management and those records maintained by the individual until re-certification is required. DO NOT SEND LESSON OR TEST TO CBSPD. For additional information regarding certification contact CBSPD - 148 Main Street, Suite C-1, Lebanon, NJ 08833 • [www.sterileprocessing.org](http://www.sterileprocessing.org).

IAHCSMM (International Association of Healthcare Central Service Materiel Management) has pre-approved this in-service for 1.0 Continuing Education Credits for a period of three years, until March 3, 2023. The approval number for this lesson is **3M-HPN 200303**.

For more information, direct any questions to *Healthcare Purchasing News* (941) 927-9345, ext. 202.

**LEARNING OBJECTIVES**

1. Discuss how vaporized hydrogen peroxide sterilization works
2. Review the keys to successful vaporized hydrogen peroxide sterilization cycles
3. Explain the recommended quality control plan for vaporized hydrogen peroxide sterilization processes

Sponsored by:



SELF-STUDY SERIES

Keys to success with vaporized hydrogen peroxide sterilization

by Craig Wallace


Low temperature sterilization processes such as ethylene oxide and vaporized hydrogen peroxide play a critical role in sterile processing departments across the United States. More and more reusable medical devices are made of advanced materials and components that cannot withstand the high temperature and moisture present in steam sterilization processes. Vaporized hydrogen peroxide (VH2O2) is the most common low temperature sterilization process used today. While this technology is not new (the first VH2O2 sterilizers appeared in U.S. hospitals in 1993), it is very complex. Understanding this process and all the requirements for successful VH2O2 sterilization are critical for patient safety.

Hydrogen peroxide (H2O2) is a very effective biocidal agent, meaning that it kills all types of microorganisms when it comes into direct contact with them. Hydrogen peroxide is a strong oxidizing agent and kills microorganisms through this mechanism. The hydrogen peroxide molecule is by nature unstable, in that it will easily “fall apart” or dissociate into

sub-components under ambient conditions. In addition, hydrogen peroxide sterilization processes use a vaporized form of hydrogen peroxide. A vapor is not the same as a true gas, in that a vapor will naturally condense on surfaces, like water vapor in the air condensing on the outside of a cold glass. Gasses, like ethylene oxide, do not condense. They are always in gas form in the temperatures and pressures found in sterilization processes. These points about hydrogen peroxide’s oxidative properties, its molecular stability, and its vapor properties are extremely important attributes that have a significant impact on how VH2O2 sterilization processes work, and why certain procedures must be followed in the sterile processing department to assure the VH2O2 processes are safe and effective.

There are many different vaporized hydrogen peroxide sterilization processes used in the various commercial sterilizers available today, but all these processes follow a similar pattern. First, there is typically a very deep vacuum step to remove air and moisture from the chamber that

Figure 1 Example of a four-pulse vaporized hydrogen peroxide sterilization process



\*STERIS V-PRO Lumen cycle1

could interfere with the process. Next, hydrogen peroxide vapor is injected into the chamber, where it will diffuse throughout a hold time. An air flush follows this hold time, then another deep vacuum is drawn to remove residual air and hydrogen peroxide. This process is considered one “pulse.” These sterilization processes will use 2-4 pulses, depending on the sterilizer and the instruments to be processed in that cycle. At the end of the process the load is aerated (washed with air) to remove residual hydrogen peroxide, or, in some systems, the load is exposed to a gas plasma process to remove the residuals. An example of a four-pulse VH<sub>2</sub>O<sub>2</sub> process is provided in Figure 1.

### Factors that can affect the VH<sub>2</sub>O<sub>2</sub> sterilization process

The chemical and physical properties of hydrogen peroxide discussed earlier (strong oxidizer, unstable molecule, vapor instead of gas) can lead to situations where a VH<sub>2</sub>O<sub>2</sub> sterilization process is not as effective as expected. Let’s take a closer look at how these properties could adversely affect a sterilization process.

1. VH<sub>2</sub>O<sub>2</sub> sterilization processes inject a fixed amount of H<sub>2</sub>O<sub>2</sub> vapor for each cycle type and for every load placed in the chamber. There are no make-ups or additions of sterilant during the sterilant exposure phase. Therefore, small loads and very large loads are exposed to the same amount of sterilant during the exposure phase. Load variation, plus the fact that the hydrogen peroxide molecule is relatively unstable and readily depletes during the exposure phase via several different chemical mechanisms<sup>2,3,4</sup> could lead to unexpected variation in the lethality of any given cycle.
2. Temperature is a critical process parameter for VH<sub>2</sub>O<sub>2</sub> sterilization. This includes the temperature of the sterilizer chamber, the medical devices, packaging, and containers when they are placed into the chamber. If the load is too cool there can be excessive condensation of the fixed amount of hydrogen peroxide vapor.<sup>5,6</sup> This can have a negative impact on the process.
3. VH<sub>2</sub>O<sub>2</sub> sterilizers are cleared by the U.S. FDA with a maximum weight limit for individual loads for each cycle type. The testing and validation supporting these clearances is completed with an understanding of the instability of the H<sub>2</sub>O<sub>2</sub> molecule, and the H<sub>2</sub>O<sub>2</sub> vapor behavior. Therefore, exceeding the weight limit for

the load can result in cycles that are not effective.

4. “Materials compatibility” typically refers to the effects of a sterilant on the medical devices. This is important for VH<sub>2</sub>O<sub>2</sub> sterilization, as hydrogen peroxide is an oxidizing chemical that can damage some materials. However, the phrase “materials compatibility” has an additional meaning when referring to VH<sub>2</sub>O<sub>2</sub> processes. In this case, it can mean the materials’ effect on the effectiveness of the process, because of interactions between the VH<sub>2</sub>O<sub>2</sub> and the materials. We have noted that VH<sub>2</sub>O<sub>2</sub> is an unstable molecule, which tends to break down into other compounds. Some materials can break down VH<sub>2</sub>O<sub>2</sub> more quickly, resulting in a loss of process effectiveness. The user must be aware that some materials (e.g. some plastics and some metals) can have a dramatic effect on the available VH<sub>2</sub>O<sub>2</sub> by absorbing, adsorbing or decomposing VH<sub>2</sub>O<sub>2</sub> at a higher rate<sup>2,3,4,7</sup>.
5. The use of extra (nonessential) materials in VH<sub>2</sub>O<sub>2</sub> sterilization is another variable that is dependent on the user and can introduce significant variation to the VH<sub>2</sub>O<sub>2</sub> sterilization process. For example, foam tray liners, polyethylene sheet tray liners, underneath guard liners, bubble wrap tray liners and tray protectors, rubber corner protectors, foam pocketed instrument protectors, CI indicator holders, transport trays, oversized disposable sterilization wrap, 600 and 650 weight disposable sterilization wrap, and preformed disposable wraps are all examples of extraneous or nonessential materials in use in healthcare facilities. As described above, because VH<sub>2</sub>O<sub>2</sub> cycles use a fixed amount of sterilant, best practices would be to limit or eliminate the use of any extra materials that could absorb the fixed amount of available VH<sub>2</sub>O<sub>2</sub> sterilant.

### The Keys to Success

These factors - a fixed amount of sterilant, a relatively unstable molecule, temperature variations, materials compatibility and use of extra nonessential materials make VH<sub>2</sub>O<sub>2</sub> processes “technique sensitive”, that is, variability introduced by the sterile processing team regarding the composition, weight, and temperature of the load can dramatically affect the outcome of the VH<sub>2</sub>O<sub>2</sub> sterilization process. In addition, there are many VH<sub>2</sub>O<sub>2</sub> sterilizer cycles available in the U.S. that have different

indications for use, different loading weight limits, different VH<sub>2</sub>O<sub>2</sub> concentrations (mg/L), and different sterilant exposure times. Let’s look at some key best practices that can help address this complicated situation by reducing variability and the potential for VH<sub>2</sub>O<sub>2</sub> sterilization process failures.

#### Chamber Loading – Weight and Spacing

Good sterilizer chamber loading practices are critical for effective VH<sub>2</sub>O<sub>2</sub> sterilization. Do not overload the chamber. Know the weight limit for your sterilizer and the sterilization cycle(s) programmed on your sterilizer. VH<sub>2</sub>O<sub>2</sub> sterilizers and cycles are cleared by the FDA with a weight limit per cycle (with exception of the STER-RAD 100S where the load weight limit is not defined)<sup>8</sup>. A good rule of thumb is to ensure there is a minimum of a hand’s width space between packages and items in the chamber (estimate about 1” space between). Items should not be stacked, should lay flat on shelves, and not contact the chamber walls or electrodes (if applicable). This type of loading allows the H<sub>2</sub>O<sub>2</sub> vapor to easily access all the packages and devices; resist the practice of adding “just one more item”.

#### Chamber Loading – Materials Compatibility

It is good practice to “know what you are loading, and only load what you know”. It is imperative that all operators of VH<sub>2</sub>O<sub>2</sub> sterilizers understand the composition of each load (devices and packaging) they place in the sterilization chamber. Some basic questions for this best practice include:

- Are the devices labeled for their specific VH<sub>2</sub>O<sub>2</sub> sterilizer model and cycle type?
- Is the total load weight below the validated and FDA-cleared weight limit?
- Is the packaging type acceptable for use in VH<sub>2</sub>O<sub>2</sub> and is the device weight under the limit for the packaging type?
- Could the device be labeled for another sterilization method like steam?
- Are there any nonessential extraneous packaging items that could be avoided?
- What is the total material composition of the load? Is the load overly weighted with items that have a higher propensity to deplete the fixed amount of VH<sub>2</sub>O<sub>2</sub>?

Understanding these basic variables for each VH<sub>2</sub>O<sub>2</sub> load will help the sterile processing team discern the effects these factors have on the process and will ultimately help assure consistent and successful process outcomes.

## Rigid Containers, Plastic Trays and Lids

Rigid containers are a very important element of VH2O2 sterilization. These containers have advantages and limitations. Their use can increase standardization of procedures and reduce waste, but they can also introduce unexpected variation over time in a VH2O2 process. Some facilities have found that rigid container surfaces and materials are designed differently e.g., some are anodized, some are not anodized, and some are expected to change in appearance over time when used in VH2O2. One manufacturer also warns against the use of soft water for the final rinse because subsequent processing in VH2O2 can cause corrosion. Containers that are not validated for VH2O2 or have worn surfaces can cause material compatibility issues with VH2O2 processes.

Plastic trays and lids are containment devices that require a sterilization wrap or pouch to maintain sterile integrity once the containment device and its contents are sterilized. Manufacturers of containment devices (including plastic trays and lids) and packaging/disposable wraps are responsible for validating that their products are compatible with VH2O2 sterilization. Plastic trays and lids must have written instructions for use and the manufacturer should provide validation data on request for the labeled sterilization modality. The Instructions for Use (IFUs) should contain the recommended maximum weight and load distribution of the containment device and its contents.

## Standards and IFUs

Vaporized hydrogen peroxide sterilization processes are complex and technique-sensitive. Consequently, understanding and carefully following all IFUs from the sterilizer manufacturer, medical device manufacturer, and packaging manufacturer are critical to success. In addition, the Association for the Advancement of Medical Instrumentation (AAMI) provides additional guidance on VH2O2 processes in AAMI ST58, "Chemical sterilization and high-level disinfection in health care facilities".<sup>9</sup> This standard provides the following points to consider for the effective use of VH2O2 sterilization:

- Follow device and sterilizer manufacturers' written IFUs
- No cellulose-based products (towels, gauze, or paper)
- Lumen sizes cleared by the FDA (based on model + cycle)
- Devices should be thoroughly cleaned and dried

- Hinged instruments should be opened
- Use only trays and mats per IFU and cleared by the FDA
- Ensure adequate sterilant contact, follow all loading recommendations
- Chemical indicators and *Geobacillus stearothermophilis* biological indicators cleared by the FDA to monitor VH2O2 sterilizers

## Quality Control

A robust quality control program is essential for VH2O2 processes. AAMI ST58<sup>9</sup> recommends the use of physical monitoring (cycle printout), chemical indicators, and biological indicators inside of Process Challenge Devices (PCDs), with the information from all three monitors combined to make a decision on the quality of the process and whether or not the devices are safe and ready for patient use. The standard provides the following guidance on the use of chemical indicators, and biological indicators in PCDs:

"A CI should be used on the outside of each package unless the internal indicator is visible"<sup>9</sup>.

"An internal CI should be used inside each package, tray, containment device (rigid sterilization container system, instrument case, cassette, or organizing tray) to be sterilized"<sup>9</sup>.

"A PCD with the appropriate BI should also be used at least daily, but preferably in every sterilization cycle (see 9.5.4.5). Each load containing implantable devices should be monitored and, whenever possible, quarantined until the results of the BI testing are available"<sup>9</sup>

The standard acknowledges that achieving the highest level of assurance for VH2O2 sterilization processes will require both consistent adherence to IFUs as well as use of a rigorous QC program:

"Most temperature sensors indicate temperature in the chamber, not at the center of packs. Improper load configuration or package composition can interfere with air evacuation and sterilant penetration, conditions that will not be revealed in the temperature recording. Therefore, physical monitoring and other indicators of sterilizer performance should never be considered a substitute for careful adherence to prescribed packaging and loading procedures."<sup>9</sup>

## Final notes regarding biological indicators for VH2O2

Currently, there is no international standard that provides performance requirements for biological indicators for VH2O2

sterilization processes. VH2O2 biological indicators from different manufacturers may be designed and tested differently, and there may be variation in the performance of these biological indicators. Since an international standard does not yet exist, the global healthcare industry has no standardization on performance requirements for BIs used in VH2O2. In the U.S., the FDA regulates biological indicators used in healthcare facilities and has a set of testing requirements for the clearance of VH2O2 biological indicators in the U.S. market. The FDA is the highest authority in the U.S. (not the sterilizer manufacturer) on the final decision on which biological indicators are cleared as compatible (safe and effective) for use in vaporized hydrogen peroxide sterilizers for healthcare facilities.

In addition, the evolution of biological indicator technology has reduced incubation times on biological indicators used to monitor VH2O2 process to 30 minutes or less. These fast readout times can enable more frequent monitoring of VH2O2 processes without interfering with the timing of the flow of instruments through the sterile processing department.

## Conclusion

Vaporized hydrogen peroxide sterilization is a complex and technique-sensitive process that requires strict adherence to IFUs from the sterilizer manufacturer, the instrument manufacturer, and the manufacturer of the packaging used in the load. In addition, a rigorous quality control program based on the use of physical sensors, chemical indicators, and biological indicators inside of PCDs is required. **HPN**

## 3M Health Care sponsored this article.

*Craig Wallace, President of Wallace Sterilization Consulting, LLC, has over 26 years of experience in the field of medical device disinfection and sterilization. Craig is the Convenor of the ISO Biological Indicator Working Group (TC 198, Working Group 4), the ISO committee responsible for international biological indicator performance standards, as well as a U.S. Technical Expert for Chemical Indicators (ISO WG 6) and Moist Heat Sterilization (WG 3). He is also the Co-Chair of the United States (AAMI) Biological Indicator Working Group, and an active member of several other AAMI working groups including chemical indicators, vaporized hydrogen peroxide sterilization, and ethylene oxide sterilization.*



References:

1. Technical Data Monograph STERIS V-PRO maX Document #10037127 Rev. C (M3635EN)
2. Block Seymour, Peroxygen Compounds. *Disinfection, Sterilization, and Preservation*, Seymour Block, 5th edition, 2001, Chapter 9.
3. U.S Patent 6,528,016 B1 March 4th, 2003 Kohler et al.
4. U.S Patent 6,875,399 B2 April 5th, 2005 McVey
5. FDA 510(k) K042116 \_ STERRAD NX Sterilizer
6. FDA 510(k) K071385\_ STERRAD 100NX Sterilizer
7. U.S Patent 5,830,409 November 3rd, 1998, Childers et al.
8. FDA 510(k) K111377 \_ Summary ASP STERRAD 100NX STERILIZER DUO CYCLE
9. ANSI/AAMI ST58, *Chemical sterilization and high-level disinfection in health care facilities*. Association for the Advancement of Medical Instrumentation. 2013 (R2018).

**CONTINUING EDUCATION TEST • APRIL 2020**

# Keys to success with vaporized hydrogen peroxide sterilization

Circle the one correct answer:

1. "Technique sensitive" is a term to describe potential variability introduced by the sterile processing department personnel that can have a significant impact on the outcome of the VH2O2 sterilization process?  
A. True      B. False
2. VH2O2 sterilization cycles and containers have loading weight limits?  
A. True      B. False
3. Best practice is to limit the use of extra or nonessential materials in VH2O2 sterilization?  
A. True      B. False
4. AAMI ST58 states: "A PCD with the appropriate BI should also be used at least daily, but preferably in every sterilization cycle"?  
A. True      B. False
5. There is no requirement for a sterilizer manufacturer to validate nor endorse indicators designed to monitor their sterilizers?  
A. True      B. False
6. Sterilizer manufacturers determine which biological indicators can be used to monitor their sterilizers.  
A. True      B. False
7. AAMI ST58 provides recommended practices for VH2O2 sterilization.  
A. True      B. False
8. All rigid containers have been validated to contain the same weight of instrumentation.  
A. True      B. False
9. The FDA clears all indicators for use in monitoring VH2O2 sterilization processes.  
A. True      B. False
10. Inappropriate loading may not be detected by the physical monitors of the sterilizer (cycle printout).  
A. True      B. False

**CONTINUING EDUCATION TEST SCORING**



The approval number for this lesson is **3M-HPN 200303**.



HPN is thrilled to now offer all CEU quizzes online by scanning the QR code on the right or visiting: <https://educationhub.hpnonline.com>. **The cost to take the quiz online remains at \$10.**

Due to rising costs if you would like to mail in your completed quiz using the method below the price is now \$50 for each test taken.

**Request for Scoring**

I have enclosed the scoring fee of **\$50 for EACH test taken** – payable to **Endeavor Business Media**. We regret that no refunds can be given. (It is not necessary to submit multiple tests separately.)

**Detach exam and return to:**

Continuing Education Division  
Healthcare Purchasing News  
2477 Stickney Point Road, Suite 221B  
Sarasota, FL 34231

